Endoscopic Dissectors

(Please Read Prior to Use)

DESCRIPTION

Endostik Endoscopic dissectors are made of an absorbent material with an embedded x-ray detectable marker that is then affixed to a 5mm shaft. Endostik Endoscopic dissectors are designed for introduction and use through all appropriately-sized trocars or a trocar with converters. The device is for use in swabbing small amounts of fluid or for blunt tissue dissection of soft tissue(s) and structures.

INDICATIONS FOR USE & KNOWN BENEFITS

Minimally invasive surgical procedures where the endoscopic dissector is used, such as, laparoscopic surgical procedures, rhinoscopy, bronchoscopy, otoscopy, cystoscopy, colposcopy, hysteroscopy, Falloposcopy and arthroscopy.

Use of endoscopic dissectors for soft tissue dissections decrease the likelihood of damage to nearby tissues or organs when used in place of sharp or metal tools.

INSTRUCTIONS FOR USE

- Step 1.** Verify count and integrity of dissector sponges when package is opened.
- Step 2.** Insert device through any appropriately-sized trocar.
- Step 3.** Use under direct visualization from an endoscope or similar method.
- Step 4.** Use the device for swabbing and blunt dissection of soft tissue(s) or structures, as needed.
- Step 5.** Verify count of sponges and dissector integrity after with drawal from the trocar/patient.

PRECAUTIONS, WARNINGS & CONTRAINDICATIONS

- Endoscopic procedures should be performed only by operating room surgeons/surgical team members having adequate training and familiarity with endoscopic techniques. A thorough understanding of operating principles, risks versus benefits, and the hazards involved in utilizing endoscopic techniques is essential to avoid the potential for injury to the patient and/or user.
- Do not use if individual packaging is damaged or open. As long as the individual package is not opened or damaged and is within the expiration date shown, the product is sterile. Use proper aseptic technique in all phases of handling.
- This product is for single use only; do not re-sterilize. Re-use can result in microbial contamination causing subsequent health deterioration of patient. Devices contain absorbent materials that cannot be cleaned. Manufacturer will not be responsible for any product that has been re-sterilized and/or reused.
- Store in a clean and dry location in manner to protect from damage or contamination.
- Tissue contact of individual devices should be limited to <1 hour.
- Count all devices before and after the procedure prior to surgical closure. In the event a device cannot be located, an x-ray can be used to locate the devices. Failure to perform count and dissector-integrity verification may result in additional procedures and/or extended surgery.
- Any serious incidents incurred during the use of the device should be reported to the manufacturer and the regulatory agency of the country in which the user/patient is established.
- Prescription medical device: Intended to be sold and used to medical practitioners for use. U.S. Federal law restricts this device to sale by or on the order of a physician.
- Used Dissectors should be considered as bio waste and be disposed of in accordance with local regulations
- Sponges may be obscured from x-ray when behind bone or in morbidly obese individuals. It is recommended that at least three views, using the optimal parameters for the imaging (x-ray) equipment, at a variety of angles (e.g., 45 degrees, 22.5 degrees, and 0 degree angles) for anterior and posterior, or the appropriate plane, be taken and examined for a missing device. If there are concerns regarding visualization, consult with your local imaging expert to establish the optimal radiographic parameters (e.g., kVp, mAs) for visualization with the imaging equipment.

Phthalate Warning

The X-ray detectable material in this device contains phthalates, specifically Diisononyl phthalate (DINP). The results of certain animal experiments have shown phthalates to be potentially toxic to reproduction. Proceeding from the present state of scientific knowledge, risks for male prenatum infants cannot be excluded in the case of long-term exposure or application. Medical products containing phthalates should be used only temporarily with pregnant women, nursing mothers, babies and infants.

Manufactured by:

First Aid Bandage Company

3 State Pier Road • New London, CT 06320 - USA
Phone: (860) 443-8499 • Fax: (860) 442-8699 • www.fabco.net

Fabco® is a registered trademark of First Aid Bandage Company
Endostik® is a registered trademark of Carwild Corporation

FABCO® is a registered trademark of First Aid Bandage Company.

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