Femoral Canal Sponge

Espanja para el canal femoral

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Femoralkanal-Schwamm

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Espanja para canal femoral

Lårbenskanalsvamp

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Femoral Kanal Süngeri

Femorale kanaalspons

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Femoral Canal Sponge

INSTRUCTIONS FOR USE

DESCRIPTION

Femoral canal sponge devices are comprised of highly absorbent PVA foam which has high tensile strength and interconnected cell structure. Each sponge has a braided tape attached for location and ease of removal from surgical site.

INDICATIONS FOR USE & KNOWN BENEFITS

Femoral canal sponge devices are intended to be used by medical professionals; operating room/surgical team members and operating surgeons used to absorb fluids during hip replacement surgery.

Use of PVA Femoral Canal Sponges will quickly and effectively dry the bone tissue while not leaving behind any fibers.

INSTRUCTION FOR USE

Step 1. Before insertion

Clean and suction the femoral canal. Introduce 10 cc of saline into the femoral canal.

Step 2. Insertion

Remove the sponge from the package and insert the unstrung end into the femoral canal. The saline will start the expansion of the sponge within the femoral canal.

Step 3. Expansion

The sponge will quickly absorb the saline and begin to expand to fill the femoral canal. A slight pressure will be exerted against the canal wall.

Step 4. Before removal

Leave the sponge in place for a minimum of one minute.

Remove the sponge prior to bone cement application by gently pulling on the drawing string or by grasping the sponge with forceps.

PRECAUTIONS, WARNINGS & CONTRAINDICATIONS

- Federal (USA) law restricts this device to use or sale by or on the order of a physician.
- Sterile. This product is for Single Use Only; do not re-sterilize. Re-use can result in microbial contamination causing subsequent health deterioration of the patient. FABCO will not be responsible for any product that has been re-sterilized or re-used.
- Do not use if individual packaging is damaged or opened. As long as the individual package is not opened, damaged or beyond its expiration date, the product is sterile. Use proper aseptic techniques in all phases of handling.
- Do not cut sponges in close proximity to the surgical procedure; fragments without x-ray detectable abilities may enter the patient's wound site. Failure to remove sponge fragments from the patient may result in a foreign body reaction.
- Store in a clean and dry location in manner to protect from damage or contamination.
- Tissue contact of individual devices should be limited to <1 hour.
- Not intended as a general surgical sponge
- There are no known side effects or residual risks related to the use of this device.
- Any serious incidents incurred during the use of the device should be reported to the manufacturer and the regulatory agency of the country in which the user/patient is established.

DISPOSAL METHOD: Used devices should be identified as biowaste or infectious waste which should be incinerated for safe disposal following local or national regulatory requirements for the disposal of biohazardous waste.

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